

From: [Shannon Day](#)
To: [David Fagan](#)
Cc: [Julie Grady](#); [Maurita Henley](#); [Lisa Kerr](#)
Subject: Newsletter Preview
Date: Tuesday, 3 November 2015 11:01:35 AM
Attachments: [image001.png](#)
[image002.png](#)
[image003.png](#)
[image005.png](#)
[image006.png](#)
[image007.png](#)
[image009.png](#)
[image011.png](#)
[image014.png](#)
[image015.png](#)

NEWSLETTER, November 2015



**Market Leaders in Prescription
Assistive Technology**

w: 1300 652 382 M: 0401 015 494
sday@motionspec.com.au



Look what's *NEW* at Motion Specialties

**Not Just a Wheelchair. A Member of the Team. Designed for Our Next Generation.
Welcome to the Most Intuitive Wheelchair You Have Seen.**

Now available

[See More](#)



***INTRODUCING* The Little Wave "Flip"**


Paediatric folding frame chair weighing 9.5kg (without rear wheels)

Designed to be as unique as your child. Flip is designed for the smallest of riders with some of the biggest advancements in the industry! From its flat folding design to its use of TAPER LOK Technology, the Flip is simply the most advanced paediatric tilt-in-space available today.



Put a flip to your fold. The Little Wave Flip integrates design elements from strollers into the durability of a wheelchair. Just because it is a wheelchair, does not mean it needs to look like one. From its flat folding design to its use of TAPER LOK technology, the Flip is simply the most advanced paediatric tilt in space available today.

[See more](#)



What's *Hot* !!

Ride Designs Java Cushion



Ride Designs' patented Java® Cushion contour works by redistributing the forces of pressure and shear from high risk areas while selectively applying the forces to areas more tolerant. The result is the potential for: • Greater sitting stability for improved functional performance • Improved skin outcomes

- Lower tissue distortion
- Lower long-term skin risk

[see more](#)

The Catalyst Difference



choosing a Catalyst.

Catalyst frames are constructed from 7000 series aluminium. A stronger and lighter vibration dampening alloy that contributes to the Catalysts superior durability and high quality construction. Whether you are in the market for a Basic or Active folding frame wheelchair you can always count on quality construction when

[Click here](#)



Did you know?

All Catalyst wheelchairs have the option to script up ready to fit solid seat bases. Special cross braces and seat hook supports free up space on the side frame for easy seating system adaptability.

[Click here](#)



For all your prescription needs contact Shannon today for a trial and assessment.

sday@motionspec.com.au

Ph: 0401 015 494

[Check It Out!](#)

Motion Specialties

3/10 Tapnor Cres

Brendale, QLD

Australia

www.motionspec.com.au

sday@motionspec.com.au



Kind Regards,

Shannon Day

Mobility Specialist

Check out what's new at Ki Mobility today!!

<http://www.kimobility.com/brochures.action>

Motion Specialties

Paragon Mobility



Ph: 0401015494

Ph: 1300 652 382



www.motionspec.com.au
www.paragonmobility.com.au



The New Little Wave CLIK Check it out Now!!!!!!

<http://www2.kimobility.com/Brochures/USA/LittleWaveCliK.pdf>

From:

Sent: Tuesday, 3 November 2015 11:00 AM

Subject: